**REF YOLA**

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***BRUGSANVISNING***

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Product code: YOLA

### 1. Characteristics

4-wheeled rollator is made of aluminium covered in paint. It has 8" diameter solid rear and front wheels, height-adjustable handles finished with anatomical grips, brakes with parking clamp, seat with height adjustable backrest and a bag. The structure is based on a crossover frame and an oval cross-section of tubes, thanks to which the resistance to axial pressure is significantly increased. The rollator is equipped with cane holder for convenient everyday use.

**The set includes:** frame with seat, front and rear solid wheels, backrest, bag, cane holder.

#### Elements of the rollator:

- A. Handgrips
- B. Parking brake
- C. Height adjustable backrest
- D. Seat with handle
- E. Bag
- F. Height adjustment
- G. Cane holder
- H. Folding mechanism for rear wheels
- I. Folding switch for front wheels
- J. Rear and front 8" solid wheels

### 2. Application

The rollator is designed for patients with reduced mobility who need stable support and help while walking. It is also recommended to the elderly as an aid in keeping balance in order to improve walking comfort and ensuring safety from falling. The use of the walker has a positive effect on the return to activity of people with locomotion and tilting to erect position problems. It guarantees better stabilization, support and security while walking.

### **3. Assembling**

#### **MOUNTING**

**Step 1.** Take the rollator and its parts out of the box.

**Step 2.** Rearrange the front wheels until you hear a characteristic “click” sound of the red folding switch (I). Pull the red folding mechanism (H) towards you and lower the rear wheels until you hear a characteristic “click” sound.

**Step 3.** Insert the handles into the openings on the frame and secure them with the locking pins.

**Step 4.** To fully unfold the rollator, grab it by the seat and pull outward. Press the seat on both sides until it clicks into place.

**Step 5.** Attach the backrest by sliding it into the special mount located under the handles. Attach the bag to the front of the rollator using straps.

**Step 6.** Attach the cane holder to the frame.

**Step 7.** Check that the rollator is assembled properly, fully unfolded, and that no parts or screws are missing. If you notice any damage or missing parts, contact the seller.

**Rollator is ready to use.**

#### **FOLDING**

- To fold the rollator, repeat above steps in the reverse order.
- Press folding switch (I) on front wheels to fold wheels upwards.

#### **BACKREST HEIGHT ADJUSTMENT**

To adjust backrest height, unscrew the cap located under the parking brake and move the backrest to desired height. Tighten the cap afterwards.

**Note!** Make sure both sides of the backrest are on the same height.

#### **HANDGRIPS HEIGHT ADJUSTMENT**

To adjust the height of the handgrip:

- Press the brake handle down to activate the parking brake and immobilize the rollator.
- Press the locking pin located on the frame and adjust height of the handle. Do the same on the other side.
- Unlock the parking brake afterwards.

**NOTE!** Both handgrips should be set at the same height.

#### **BRAKE FUNCTION**

- To stop the rollator, pull the brake handle up, then hold.
- When not in motion or when using the seat, lock the parking brake with the parking clamp by pulling the brake handle down.

- To release the parking brake, pull the handle upward.

#### 4. Technical data

Overall width	65 cm
Overall length	69 cm
Seat width	46 x 25 cm
Seat height from the ground	55 cm
Handgrip height adjustment	81-97 cm
Rollator weight	8,20 kg
Maximum capacity	136 kg
Maximum capacity of the bag	5 kg
Frame colour	matte black, turquoise, grey

#### 5. Notes

**Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.**

#### 6. Maintenance

- Clean the frame with damp cloth and dry thoroughly.
- Do not use strong polishing agents due to the damage of protecting coatings.
- Clean the dirt from moving parts of the rollator regularly (hinges, clamping screws, locks).

#### 7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

#### 8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

#### 9. Labels

	Reference number		Manufacturer
---	------------------	---	--------------

	LOT number		UDI code
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: [www.timago.com](http://www.timago.com). Thank you for choosing Timago!

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