FLASH-TIM

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USER MANUAL

Electric wheelchair

Product code: FLASH-TIM

1. Security information

- Read the manual before first use.
- Do not use the wheelchair until fully assembled.
- Use of the wheelchair by people who are unable to operate it independently is not recommended.
- Maximum load capacity must not exceed 120 kg.
- The wheelchair can only be used by one person at a time.
- It is recommended to keep your feet on the footrest and your forearms on the armrests when driving.
- Always check the battery level on the control panel before use.
- This electric wheelchair should only be used on a flat surface or at low elevations. It should not be ridden on a surface inclined at an angle of more than 8 degrees.
- Do not drive the wheelchair at maximum speed when going uphill or downhill.
- Do not use the wheelchair on soft and slippery surfaces such as mud, sand etc.
- Do not ride the wheelchair sideways on highly elevated areas.
- Do not put any additional load such as a trailer etc. on the wheelchair.
- Do not repair or modify the construction of the wheelchair by yourself. Use only original spare parts provided by the manufacturer.
- Do not stand up or sit down on the wheelchair when the control module is switched on.
- Do not use the wheelchair without anti-tip wheels.
- When the freewheel (manual) mode is on, do not get on or off the wheelchair without assistance.
- Do not lean forward to grab items laying in front, on the side or behind the wheelchair. Do not stand on the footrest. It might result in the wheelchair falling over.
- Avoid rapid changes of the direction while driving on elevations.
- Always slow down while driving downhill or making turns.
- Before use check if the wheels are correctly mounted.
- Do not make sudden moves while using the joystick.

- When storing, protect the wheelchair from rain and adverse weather conditions.
- Before turning on the wheelchair, make sure that the levers are set to electric mode.
- Avoid areas with horizontal gaps.
- If the wheelchair is not going to be used for a long period of time – remove the battery.
- Do not use the wheelchair on the road. It might result in a car accident.
- Do not store the wheelchair in places with high temperatures and humidity.

2. Product description and application

Electric wheelchair is designed for people who have problems with mobility, helping disabled patients with their everyday tasks. It can replace a standard wheelchair.

Electric wheelchair is suitable for patients with:

- disability, who have problems with unassisted movement,
- dysfunctions and reduced movement,
- heart and cardiovascular system diseases,
- weakening body due to advanced age.

Medical device incidents

Any serious incident involving a medical device must be reported to the manufacturer and to the competent authority of the European Union member state where the user or patient resides.

Icon description

Symbol	Meaning
	Safety Classification type B
	Important! Read the user manual
IP44	Protection against access to hazardous parts with a tool, protection against foreign solids with a diameter of 1 mm and larger, protection against water spray from any direction
	Fragile
	This side up

Protect from humidity and moisture

A. ELEMENTS OF THE WHEELCHAIR

1. Handle
2. Backrest
3. Seat
4. Armrest
5. Control module
6. Support wheel for transportation
7. Rear wheel 12"
8. Front wheel 8"
9. Footrest

10. Bag
11. Security buckle
12. Mounting screw
13. Folding button
14. Lithium battery
15. Engine (mode) levers

B. PRODUCT CHARACTERISTICS

- Aluminium, folding frame
- Solid rear 12" and front 8" wheels
- Control module: on/off button, steering joystick, horn button, speed adjustment, battery level indicator, mute button
- Safe and reliable electromagnetic brake
- Fully closed, maintenance-free battery
- Security belt and anti-tip wheels
- Pocket for small goods located at the back
- Two modes of operation: electric and manual
- Strap - for securing a folded wheelchair

C. TECHNICAL PARAMETERS

1. **Product type:** for indoor and outdoor use
2. **Maximum speed:** ≤ 6 [km/h]
3. **Maximum capacity:** 120 [kg]
4. **Braking distance on horizontal surface:** $\leq 0,8$ [m]
5. **Max. range:** 20-25 [km]
6. **Battery:** lithium, 24 [V] x 12 [Ah]
7. **Max. angle of elevation:** $\leq 10^\circ$
8. **Min. turning radius:** 1,20 [m]
9. **Wheelchair weight:** 19,10 kg
10. **Battery charging time:** 6 hours
11. **Battery weight:** 3,30 kg

Note! Some parameters might vary depending on user weight, road conditions and battery usage.

- **Safety Classification:** type B
- **Ingress protection level for engine:** IP44
- **Engine power level:** 2x150[W]
- Wheelchair cannot be used in areas with flammable gas mixed with air, oxygen or nitrous oxide.

3. Assembly and adjustment

- **Unfolding:** Take out electric wheelchair from a box. Hold the seat frame with one hand and backrest frame with another, spreading them apart till fully unfolded (characteristic "click" sound; safety buckle should be locked). Pull down the footrest.

- **Folding:** To fold the wheelchair, press on the safety buckle located behind the wheelchair (away from you) and slowly fold wheelchair by moving backrest forward. Fold footrest.
Note! Disconnect the power cable from the battery before unfolding.
- **Regulation of control module position:** Regulation should be only done while sitting in wheelchair. Loosen mounting screw and move mounting bracket forward or backward to set it to appropriate length, then tighten the screw.
- **Armrest tilting:** Press the buttons located on the outside of the wheelchair to lift the armrests. To lower the armrest, pull them down until you hear a “click” sound.
- **Anti-tip wheels assembly:** Put the wheels into the mounting openings located at the back of wheelchair and pull the lever next to it to secure the wheels in place.
- **Removing the battery:** to remove the battery, unlock the buckle located underneath and slowly pull the battery out.

4. Operation

A. CONTROL MODULE FUNCTIONS



- A. Joystick
- B. Power on/off button
- C. Speed level indicator
- D. Horn
- E. Speed adjustment buttons
- F. Battery indicator
- G. Mute button
- H. USB port and module with flashlight

- **Joystick (A)** – used to control speed and direction (forward, backward and turn). Ride direction depends on director in which joystick is pushed. The

further joystick is tilted from middle of control panel the faster wheelchair moves. Wheelchair will automatically lose speed when joystick will be released. During wheelchair stop do not move joystick rapidly.

- **Speed indicator has 5 levels (C):** 1 – slowest, 5 – fastest.
- **Mute button (F)** – mutes the assistant voice (only available in English).
- **Horn (D)** – press joystick to use horn.

LED display	Status	Meaning/Description
	All LED lights are turned off	Wheelchair is turned off or in sleep mode.
<p>1.</p>  <p>2.</p>    <p>3.</p> 	<p>1. All lights are on</p> <p>2. Few lights are on</p> <p>3. Only red light is on</p>	<p>When the wheelchair is turned on, the number of lights shows the battery charge level:</p> <p>5 lights = 81-100%</p> <p>4 lights = 61-80%</p> <p>3 lights = 41-60%</p> <p>2 lights = 21-40%</p> <p>1 light = ≤ 20%</p> <p>If there's only one light (red) left, battery is running low. Charge battery as fast as possible.</p>

B. MANUAL AND ELECTRIC MODE SWITCHING

To use the wheelchair in manual mode, move down the left and right engine levers located at the back. To use it in electric mode, move the levers up.

Electric mode

Manual mode

- Set includes a charger dedicated to the battery with power of 24 [V].

- Plug the charger into the charging port located in the control module or on the battery.

NOTE: Before plugging the charger in, turn off the control panel. If connecting charger to battery directly, unplug battery cables before charging.

- Orange light indicates that the battery is being charged, while green light informs that it is fully charged. After the charging is finished, the charger turns off to avoid overcharging.
- Check if circuit connection is correct and not damaged.
- **Do not** modify charging circuit connection.
- To avoid fire or burnout, do not unplug charging circuit while battery is charging.
- To check the current battery level on control module, plug off the charger.

D. DIRECTIONS FOR USE

- Select the desired driving mode.

Note! Do not switch the mode while driving on elevation.

- Push the power button to check the proper functioning of brakes. Brakes are working correctly (electric mode) if you cannot push the wheelchair. Otherwise, contact the Seller.
- Tilt the footrest and sit comfortably on the seat. Lower the footrest back to horizontal position.

Note! Do not use the footrest as a support when sitting down or standing up from the wheelchair. This may cause the wheelchair to tilt and fall.

- Sit on the wheelchair and turn the control panel on. The LCD panel will light up. Joystick should be in neutral position.
- Take caution when using for the first time. Control the speed and direction carefully. The electromagnetic brake will get released only after pushing the joystick (electric mode). To speed up, push the joystick further away from its neutral position. To slow down – do the opposite.
- To stop the wheelchair, release the joystick. Pressing the power button will also stop the wheelchair.
- Using dedicated buttons on the control panel, adjust the wheelchair's speed according to user's physical state and road conditions.
- Wheelchair should only be used on flat surfaces. Driving on muddy or uneven surfaces can cause control system damage.

5. Technical data

Product code	FLASH-TIM
Overall height	93 cm
Overall width	62 cm
Overall length	90 cm
Seat height from the ground	54 cm (with cushion)
Backrest height	41 cm
Seat cushion dimensions	45 x 45 x 6 cm
Battery weight	3,30 kg
Maximum capacity	120 kg
Wheelchair weight*	19,10 kg*
Wheels	front: solid, 8" rear: solid, 12"
Folded width	41 cm**

* wheelchair weight without battery

** measured without cushion attached

6. Cleaning and maintenance

- It is recommended to regularly inspect the construction and its elements (wheels, frames, cables, control module, brakes etc.).
- Any used or damaged wheelchair parts need to be replaced (in such case, please contact the Seller).

7. Battery maintenance and cleaning:

- 1) Regularly check the battery level displayed. If the level is very low, charge it immediately to avoid a voltage drop (which may affect battery life).
- 2) If wheelchair won't be used for longer period of time, charge battery at least once a month.
- 3) Battery has a certain lifespan. Exchange your battery if it wears out.
- 4) Battery should be kept clean and away from kids. Avoid causing any shock or damage to the battery.
- 5) Charge it regularly to expand its lifespan.
- 6) Used batteries should be utilized according to local environmental regulations.

8. Storage and transport conditions

A. TRANSPORTATION

- Take caution during transportation. Do not throw, turn over or put pressure on the product and its packaging.
- For wheelchair folding, check „Assembly and adjustment”.

B. STORAGE

Wheelchair should be kept in dry and well-ventilated space. Avoid: high temperature and its rapid changes, acids, alkalis and corrosive chemicals.

C. STORING CONDITIONS

Temperature: -25[°C] ~+55[°C]

Humidity: 30-75%

9. Troubleshooting

- In the event of a malfunction, the control module will indicate it with a beep and a flashing diode. The number of beeps and LED flashes will help to identify the problem.

Number of flashes	Type of problem	Cause / possible solution
1	User error	User exceeded the limit set by controller. Wheelchair does not turn on.
		Release the joystick back to center position and try again.
2	Interface error	Control panel interface error occurs.
		Contact with the Seller to replace the panel.
3	Battery/charger failure	Charger or battery may have been damaged.
		Check if the power cords are tightened.
4	Left engine failure	Check if power cords connected to the left engine are tightened.
5	Right engine failure	Check if power cords connected to the right engine are tightened.
6	Failure of left electric brake	Check if power cords connected to the left brake are tightened. Look for damages or connection failures.

7	Failure of right electric brake	Check if power cords connected to the right brake are tightened. Look for damages or connection failures.
8/9	Control module error	Check the connection. Please contact the Seller.

10. Warranty information

- All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained. As user-friendly company, we provide only tested products, in terms of material, quality and functionality. If you have any other question regarding service procedures, contact us.
- **Warranty does not covers:**
 - **wearable parts:** upholstery, tire, armrests, battery
 - **damages resulted from:** self-repairs, changes in product construction, drop, transport, inappropriate usage, noncompliance to the instruction, natural disasters (fire, earth shake, floods)
- Product code is located on wheelchair frame label.

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Anhui Gmed Technology Co., LTD.
SUNGO Europe B.V.